

# GLUTARALDEHYDE

## WARNING!

**CAN BE IRRITATING AND HARMFUL TO SKIN,  
EYES, OR IF INHALED**

**May cause allergic skin reaction, cough, headache,  
severe eye irritation.**

**PRECAUTIONARY STATEMENTS:** Avoid breathing vapor. Use with adequate ventilation. Wear splash goggles and/or face shield, gloves, and protective clothing.

**FIRST AID PROCEDURES:** If inhaled- Remove to fresh air; resuscitation, oxygen, medical attention may be necessary. Eye contact - Flush with flowing water for at least 15 minutes. Seek medical attention. Skin contact - remove contaminated clothing, flush area with water for 15 minutes. Ingestion - Do not induce vomiting. GET MEDICAL ATTENTION.

**Consult MSDS for further health and safety information.**

CAS NO. [111-30-8]

Brand Name Used \_\_\_\_\_

**Kem Medical Products Corp.**

**(800) 553-0330**