

Price List

GUS SOAK STATIONS

Model #	Retail Price	Description
Ultrasound		
vaginal – rectal - transesophageal		
G10VP	\$ 789	Wall-mounted station for one vaginal and/or rectal probe. ONE soak tube, one rinse tube.
G14KA	\$ 1,130	Wall-mounted station for two vaginal and/or rectal probes. TWO soak tubes, one rinse tube.
G14T-3	\$ 1,425	Wall-mounted station for two transesophageal probe. TWO straight soak tubes, one rinse tube.
G14TC-3	\$ 1,747	Wall-mounted station for two transesophageal probes. TWO curved soak tubes, one rinse tube.
Small Scopes		
hysteroscopes – cystoscopes – ENT scopes		
G10ENT	\$ 789	For partial immersion of small scopes . ONE soak tube, one rinse tube. Wall mounted only.
G10-U2	\$ 852	For partial immersion of one Cystoscope/Hystocope . One soak tube, one rinse tube. Wall mounted only.
G14ENT-3	\$ 1,130	For partial immersion of small scopes . TWO soak tubes, one rinse tube. Wall mounted only.
G14-U3	\$ 1,306	For partial immersion of two Cystoscopes/Hystoscopes . Two soak tubes, one rinse tube. Wall mounted only.
G14FH	\$ 1,130	General purpose hood . For small instruments and accessories. (14.5" wide X 16.5" high)
G14-U2-AC	\$ 1,192	For partial immersion of one Cystoscope/Hystoscopes and one US Probe . Two soak tubes, one rinse tube. Wall mounted only.
G17	\$ 2,136	For full immersion of small scopes (Space Saver Station* - 18" wide X 16" high)
G18-E	\$ 2,493	For full immersion of small scopes (Elevated Countertop Station - 19" wide X 24.5" high)
Large Scopes		
hysteroscopes – cystoscopes – ENT scopes		
G25-E	\$ 2,549	Countertop station (27" wide x 27 1/2" deep x 24 1/2" high)
G25-S	\$ 2,379	Space Saver Station (27" wide x 27 1/2" deep x 16 1/2" high)
G32-E	\$ 2,833	Countertop station (33 1/4" wide x 27 1/2" deep x 24 1/2" high)
G32-S	\$ 2,674	Space saver station (33 1/4" wide x 27 1/2" deep x 16 1/2" high)
G36-E	\$ 3,742	Space saver station (32" wide x 43" deep x 24 1/2" high)
G36-S	\$ 3,526	Space saver station (32" wide x 43" deep x 16 1/2" high)
G44-E	\$ 3,895	Countertop station (45" wide x 27 1/2" deep x 24 1/2" high)
G50-E	\$ 4,633	Countertop station (54" wide x 27 1/2" deep x 24 1/2" high)

S = Space Saver Stations are designed to fit under upper wall cabinets
E = Elevated – gives extra work height



75 Price Parkway, Farmingdale, NY 11735
Telephone: (631) 454-6565 • (800) 553-0330
Email: mail@kemmed.com

FAX# (631) 454-8083
Web Address: www.kemmed.com

January 1, 2011

Part #	Description	Single Filter	Price All at Once 2 Year Plan*	% Discount
F10	For G10 models	\$ 181.00	\$ 543.00	25%
F141	For G14 models	\$ 239.00	\$ 717.00	25%
F171	For G17 models	\$ 250.00	\$ 750.00	25%
F181	For G18 models	\$ 250.00	\$ 750.00	25%
F201	For G20 models	\$ 301.00	\$ 903.00	25%
F221	For G22, SSD & DSD models	\$ 272.00	\$ 816.00	25%
F251	For G25 models	\$ 296.00	\$ 888.00	25%
F321	For G32 models	\$ 312.00	\$ 936.00	25%
F361	For G36 models (dual pack)	\$ 477.00	\$ 1,431.00	25%
F181-251	For G44 models (dual pack)	\$ 522.00	\$ 1,566.00	25%
F50	For G50 models (dual pack)	\$ 545.00	\$ 1,635.00	25%

* Every six months, you will be sent a letter reminding you to change the filter

Payment terms: Net 30 or credit card. Shipping and handling charges prepaid and added to invoice. Returns require an RMA # and are subject to a 25% restocking fee plus shipping and damage costs. Prices are subject to change without notice.

Specifications for GUS™ Ductless Fume Hoods

January 1, 2011

Model	Air Treated (CFM)	FPM w/door closed low	FPM w/door open high	Voltage	Current w/door closed	Current w/door open	Dimensions (outside) H x W x D	Dimensions (inside) W x D
G10VP	42	200	85	110V 50/60Hz	18 watts	18 watts	22.25"H x 10"W x 10.25"D	Not Applicable
G10ENT	42	200	85	110V 50/60Hz	18 watts	18 watts	31.5"H x 10"W x 10.25"D	Not Applicable
G14ENT-3	50	200	N/A	110V 50/60Hz	19 watts	19 watts	16"H x 15"W x 12.25"D	Not Applicable
G14ENT-AC	50	200	N/A	110V 50/60Hz	19 watts	19 watts	16"H x 15"W x 12.25"D	Not Applicable
G14FH	50	60	N/A	110V 50/60Hz	19 watts	19 watts	16.5"H x 14.5"W x 21.25"D	14"W x 16.25"D
G14KA	50	200	75	110V 50/60Hz	19 watts	19 watts	20"H x 15"W x 12.25"D	Not Applicable
G14T3	50	N/A	N/A	110V 50/60Hz	19 watts	19 watts	48"H x 14.5"W x 15"D	Not Applicable
G14TC-3	50	N/A	N/A	110V 50/60Hz	19 watts	19 watts	48"H x 14.5"W x 15"D	Not Applicable
G17	110	120	65	110V 50/60Hz	19 watts	38 watts	16.5"H x 18.25"W x 24.5"D	17.5"W x 23.56"D
G18E	110	120	65	110V 50/60Hz	19 watts	30 watts	24"H x 19.5"W x 24.5"D	18"W x 23.56"D
G25-E	210	120	78	110V 50/60Hz	40 watts	75 watts	24.5"H x 27"W x 27.5"D	25.6"W x 22"D
G25-S	210	120	78	110V 50/60Hz	40 watts	75 watts	16.5"H x 27"W x 27.5"D	25.6"W x 22"D
G32-E	220	120	68	110V 50/60Hz	40 watts	75 watts	24.5"H x 33.25"W x 27.5"D	32"W x 22"D
G32-S	220	120	68	110V 50/60Hz	40 watts	75 watts	16.5"H x 33.25"W x 27.5"D	32"W x 22"D
G36-E	220	120	68	110V 50/60Hz	40 watts	75 watts	24.5"H x 37.75"W x 27.5"D	37"W x 22"D
G36-S	220	120	68	110V 50/60Hz	40 watts	75 watts	16.5"H x 37.75"W x 16.5"D	37"W x 22"D
G44-E	300	120	69	110V 50/60Hz	66 watts	112 watts	24.5"H x 45"W x 27.5"D	44" W x 22"D
G50-E	400	120	45	110V 50/60Hz	80 watts	150 watts	24.5"H x 54"W x 27.5"D	51"W x 22"D



75 Price Parkway, Farmingdale, NY 11735
Telephone: (631) 454-6565 • (800) 553-0330 FAX# (631) 454-8083
Email: mail@kemmed.com Web Address: www.kemmed.com

ACCESSORIES

January 1, 2011

Part #	Retail Price	Description
Soaking Cups/Tubes		
TUBE-C	\$ 173.00	Curved tube for G14TC station
TUBE-VP2	\$ 33.00	12" Replacement Container (2-pack)
TUBE-VP3	\$ 48.00	12" Replacement Container (3-pack)
TUBE-2.5	\$ 163.00	2.5" Soak Tube (for larger probes)
Trays		
TR4-6A/2	\$ 37.00	4" x 6" x 2" Tray w/lid
TR6-10C	\$ 38.00	6" x 10" x 6" non-autoclavable tray with lid
TR6-10A	\$ 60.00	6" x 10" x 6" autoclavable tray with lid
TR13-21C/6	\$ 60.00	13" x 21" x 6" non-autoclavable tray with lid
TR13-21A	\$ 178.00	13" x 21" x 6" autoclavable tray with lid
TR15-21W	\$ 125.00	15" x 2" x 6" autoclavable tray with lid
TR-HS-U	\$ 130.00	24" x 8.75" x 5.5" non-autoclavable tray with lid
Carts		
CT25	\$ 487.00	Cart for G25 stations
CT32	\$ 514.00	Cart for G32 stations
Other		
AC-10-S	\$ 60.00	Stand for G10VP
AC-14-S	\$ 70.00	Stand for G14KA
DS02	\$ 416.00	Stand-alone disposal system
AC-INLF	\$ 45.00	Five (5) in-line filters for disposal system
AC-TEST	\$ 163.00	Leak Test Assembly
AC-G14RF	\$ 375.00	Replacement front
AC-BRUSH	\$ 60.00	Tube Brush
AC-CLIP-5	\$ 27.00	Grip Clips (5-pack)
AC-SUTC	\$ 195.00	Wall-mounted storage tube for TEE Probes
AC-TCRK	\$ 319.00	Probe storage rack with bracket