

OPA (o-Phthalaldehyde)

WARNING!

**CAN BE IRRITATING AND HARMFUL TO SKIN,
EYES, OR IF INHALED**

**May cause allergic skin reaction, cough, headache,
severe eye irritation.**

PRECAUTIONARY STATEMENTS: Avoid breathing vapor. Use with adequate ventilation. Wear splash goggles and/or face shield, gloves, and protective clothing.

FIRST AID PROCEDURES: If inhaled- Remove to fresh air; resuscitation, oxygen, medical attention may be necessary. Eye contact - Flush with flowing water for at least 15 minutes. Seek medical attention. Skin contact - remove contaminated clothing, flush area with water for 15 minutes. Ingestion - Do not induce vomiting. GET MEDICAL ATTENTION.

Consult MSDS for further health and safety information.

Kem Medical Products Corp.

(800) 553-0330