



## KemSure® Spill Control Kits

### Overview:

KemSure® Spill Control Kits provide all components necessary to rapidly respond to small spills of hazardous liquid chemicals. Kits are available for use with spills of Formaldehyde, Glutaraldehyde, OPA and OPA-C.


Chemsorb™ Spill Response Pillows absorb and collect the liquid residue, providing a fast-acting spill solution. KemSafe™ Neutralizing Powders deactivate and neutralize hazardous chemicals, rendering them non-toxic for safe disposal. Personal protective equipment in kits provides for safe handling.

A KemSure Universal Spill Control Kit is also available for spills that do not require neutralization, or for which a neutralization method is not available.

### Feature Highlights:

- Available in packaging choice of pail or tote bag for storage convenience.
- Color-coded for easy identification of corresponding chemical.
- Convenient and cost-effective way to dispose of small liquid chemical spills.

### About Kem Medical Products Corp.

We offer a complete line of spill control kits, as well as a full range of health and safety products and services for healthcare facilities. Kem Medical Products Corp. has been proudly serving the healthcare industry since 1980. Our mission is to provide healthcare professionals with the high quality products and solutions they need in order to maintain a safe and effective work environment. 



### Kem Medical Products Corp.

400 Broadhollow Road, Suite 2, Farmingdale, NY 11735 ■ +1.800.553.0330 ■ [www.kemmed.com](http://www.kemmed.com)



## KemSure® Spill Control Kits: Performance and Item Details

### KIT CONTENTS

- One pair Nitrile Gloves
- One AAMI Level-3 Isolation Gown
- One Procedure Mask with Elastic Earloops
- One Full-Length Face Shield
- One Pair Non-Skid Knee-High Boot Covers
- One Lightweight Bouffant Cap
- Two Spill Response Pillows
- Neutralizing Powder (if applicable)
- Temporary Disposal Bag with Label
- One Inventory Control/Refill Card
- IFU/SDS Documents

### PRODUCT OFFERINGS

#### KemSure Spill Control Kits

Item	Packaging	Kit Type	Summary
9500-P	Pail	Universal Spill Control Kit	Each spill response pillow absorbs 1/4-gallon of liquid. Two pillows included.
9500-T	Tote		
9520-P	Pail	Formaldehyde Neutralizing Spill Control Kit	Each neutralizing pillow absorbs 1/4-gallon of liquid. Two pillows included.
9520-T	Tote		
9540-P	Pail	OPA-C Neutralizing Spill Control Kit	Powder neutralizes 5-liters of OPA-C liquid. Each spill response pillow absorbs 1/4-gallon of liquid. Two pillows included.
9540-T	Tote		
9560-P	Pail	OPA/Glutaraldehyde Neutralizing Spill Control Kit	Powder neutralizes up to 2-gallons of OPA-C liquid. Each spill response pillow absorbs 1/4-gallon of liquid. Two pillows included. <i>*Extra pillows or larger pillows may be purchased.</i>
9560-T	Tote		
9540-S	Boxed Set (8 containers per Kit)	OPA-C Residual/Surplus Neutralizing Kit (for residual/surplus liquid remaining in used OPA-C sterilant containers)	Each container (8 per kit) neutralizes 2.28 oz. of residual OPA-C liquid.

#### Additional Products

Item	Description
9520-S	Two (2) Formaldehyde Neutralizing Pillows (each for 1/4-gallon Formaldehyde liquid).
9540-R	OPA-C Kit Neutralizing Powder Refill (for up to 5-liters of OPA-C liquid).
9560-R	OPA/Glutaraldehyde Kit Neutralizing Powder Refill (for up to 2-gallons of OPA-Glutaraldehyde liquid).
9500-RP	<i>*Two (2) Spill Response Pillows (each absorbs up to 1/4-gallon of liquid).</i>
9251	<i>*Six (6) Large Spill Response Pillows (each absorbs up to 1/2-gallon of liquid).</i>